

53 88

EC

HACH®

53 88

**Key points. Hauptteile.
Puntos clave. Points clé.
Punti chiave. Pontos chaves.**

① MP-5 connector.
MP5-Stecker.
Conector MP-5.
Connecteur MP-5.
Connettore MP-5.
Conector MP-5.

② PG 13.5 thread.
Montagegewinde PG 13,5.
Rosca de sujeción PG 13,5.
Filetage PG 13,5.
Filettatura PG 13,5.
Rosca de fixação PG 13,5.

③ Titanium body, electrode 1.
Titangehäuse, Elektrode 1.
Cuerpo de titanio, electrodo 1.
Corps en titane, électrode 1.
Corpo in titanio, elettrodo 1.
Corpo em titânio, eletrodo 1.

④ Pt 1000 temperature sensor.
Temperatursensor Pt 1000.
Sensor de temperatura, Pt 1000.
Capteur de température, Pt 1000.
Sensore di temperatura, Pt 1000.
Sensor de temperatura, Pt1000.

⑤ Electrode 2, titanium.
Elektrode 2, Titan.
Electrodo 2, titanio.
Electrode 2, titane.
Elettrodo 2, titanio.
Eletrodo 2, titânio.

- * Install the cell in immersion housings
(see mounting).
*Installation in Taucharmaturen
(siehe Montage).*
Instalar en sondas de inmersión
(ver montaje).
*Pour installer en sondes d'immersion
(voire montage).*
Installare in sonda a immersione
(ved. montaggio).
*Instalar a célula no acessório para imersão
(veja a seção Montagem).*

(EU Directive 2002/96/EC)

Conductivity cell

Introduction

The **53 88** is a cell for reservoirs and tanks. They include an integrated Pt 1000 temperature sensor.

SAFETY NOTE

Before use with hazardous substances, find out about the relevant hazards, protective measures, and what to do in an emergency.

Cell preparation

Before use for the first time immerse the electrodes in ethanol for 15 s, and rinse with distilled water. Apply the same procedure for a cell not being in use for long time.

Installation

The cell should always be mounted in immersion or flotation housings, with PG 13,5 thread.

 See *Housing instruction manual*.

Cable. Connected to multiple cables with MP-5 connector.

Calibration

Before measurement, calibrate with standard solutions.

 See *User's Manual of conductivity meter*.

Measurements. Recommendations

- Both, calibration and measurement must be performed under slight stirring of standards and samples. The stirring speed in the two cases should be identical.
- The sample must cover the measuring electrodes ⑤.
- Ensure that during measurement there are not air bubbles over the measuring electrodes.
- After measurement, rinse the cell with distilled water or with sample to be measured.

 With time and use the titanium alters its colour. This change does not affect the good functioning of the cell.

Maintenance

Keep the electrodes surface clean ⑤.

If the thorough rinse with distilled water is not sufficient, suitable solvent for the elimination of sample rests can be used. This solvent must not produce chemical attack over the cell materials.

Inorganic contamination: Immerse the cell in diluted hydrochloric acid or in chromic acid. Afterwards rinse thoroughly with distilled water.

Organic contamination: Use suitable solvent (see cell's materials), next utilize alcohol and distilled water as in Cell preparation is described.

Storage

Store the cell dry.

Life of the cell

The cell lifespan is unlimited if suitable maintenance is applied and the cell does not break.

Troubleshooting

... Reading 0 independently of the measured sample.

Broken cable.

Replace the cell.

... Unstable reading.

... Impossible calibration.

Dirty measuring electrodes or air bubbles on their surface.

Clean or shake slightly the cell.

... The obtained readings differ a lot from the real values.

Verify that the conductivity is covered by the measuring interval of the cell.

Re-calibrate using standards in good conditions.

... slow cell response

Regenerate the cell.

Specifications

Approx. constant	0,6 cm ⁻¹
Measuring range	10 µS/cm to 80 mS/cm
Operating temp.	0 to 80 °C
Max. pressure	6 bar (25 °C)
Temperature sensor	Pt 1000
Body material	titanium
Electrode material	titanium
Thread	PG 13,5
Connector	MP-5
Min. immers. depth	10 mm 

Applications

Drinking waters, municipal waters, water for irrigation, industrial waters, wastewater, etc.

CERTIFICATE OF QUALITY

The attached cell **53 88** has passed the quality test and fulfils the following specifications:

- EC error ≤ 5 %
between 10 µS/cm and 80 mS/cm.
- Temperature error ≤ 0.5 °C.

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